

HI 3811 Alkalinity Phenolphthalein and Total

Alkalinity is the quantitative capacity of a water sample to neutralize an acid to a set pH. This measurement is very important in determining the corrosive characteristics of water due primarily to hydroxide, carbonate and bicarbonate ions. Other sources of alkalinity can be from anions that can be hydrolyzed such as phosphates, silicates, borates, fluoride and salts of some organic acids. Alkalinity is critical in the treatment of drinking water, wastewater, boiler & cooling systems and soils.

Alkalinity Conversions

 $1 \text{ meq/L} = 50 \text{ mg/L CaCO}_3 = 2.8 \text{ dKH}$ $1 \text{ mg/L CaCO}_3 = 0.02 \text{ meg/L} = 0.056 \text{ dKH}$ $1 \text{ dKH} = 0.36 \text{ meg/L} = 17.86 \text{ mg/L CaCO}_3$

There are three methods of expressing alkalinity generally used:

mg/L CaCO₃ = milligrams of CaCO₃ per liter water

200

865 g

meq/L = milliequivalents per liter

dKH = degrees of carbonate hardness

0,3 0,1

METHOD	RANGE*	SMALLEST INCREMENT	CHEMICAL METHOD	# TESTS	WEIGHT	
HI 3811 Alkalinity (as CaCO ₃) Phenolphthalein and Total						
titration	0-100 mg/L (ppm) 0-300 mg/L (ppm)	1 mg/L (ppm) 3 mg/L (ppm)	phenolphthalein/ bromphenol blue	110 avg.	460 g	
UI 2004 & Allia Barta Tarad						
HI 38014 Alkalinity Total						
titration	0-500 gpg	5 gpg	bromphenol blue	100	363 g	
HI 38013 Alkalinity, Phenolphthalein and Total						
titration	0.0-10.0 gpg	0.1 gpg	phenolphthalein/	200	865 a	

0.2 gpg

bromphenol blue

titration

HI 3811 Alkalinity Test Kit

The HANNA alkalinity test kit makes monitoring easy, quick and safe. The compact size gives the user the versatility to use the kit anywhere. The design makes the kit easy to handle and, except for alkalinity titrant, practically prevents accidental injury or damage due to spills.

Alkalinity can be measured as phenolpthalein alkalinity and total alkalinity. The phenolpthalein alkalinity is determined by neutralizing the sample to a pH of 8.3 using a dilute hydrochloric acid solution, and a phenolpthalein indicator. Since bicarbonate ions can be converted to carbonic acid with additional hydrochloric acid, the phenolpthalein alkalinity measures total hydroxide ions, but only half of the bicarbonate contribution.

HI 38013 Phenolphthalein and Total **Alkalinity Test Kit**

In this case the phenolphthalein alkalinity is determined by neutralizing the sample to a pH of 8.3 using a dilute sulfuric acid solution and a phenolphthalein indicator. Since bicarbonate ions can be converted to carbonic acid with additional sulfuric acid, the phenolphthalein alkalinity measures total hydroxide ions, but only half of the carbonate contribution.

HI 38014 Total Alkalinity Test Kit

Total alkalinity is determined by neutralizing the sample to a pH of 4.5 using a dilute sulfuric acid solution and a bromophenol blue indicator.

ORDERING INFORMATION

HI 3811 test kit comes with 10 mL phenolpthalein indicator, 10 mL bromophenol blue indicator, 120 mL alkalinity titrant, 10 mL calibrated vessel, 50 mL calibrated vessel, and calibrated syringe with tip. HI 38013 test kit comes with 10 mL phenolphthalein indicator, 10 mL bromophenol blue indicator, 105 mL alkalinity reagent (2), 20 mL calibrated plastic vessel with cap and 1 mL syringe with tip. HI 38014 test kit comes with 10 mL bromophenol blue indicator, 110 mL total alkalinity reagent, 20 mL calibrated vessel with cap and 1 mL syringe with cap.

ACCECCODIEC

ACCESSORIES					
HI 3811					
HI 3811-100	Spare reagent for 100 tests				
HI 38013					
HI 38013-100	Spare reagent for 100 tests				
HI 38014					
HI 38014-100	Spare reagent for 100 tests				



0.0-20.0 gpg

^{* 1} gpg = 17 ppm CaCO;