



KYORITSU

PACKTEST
ION SELECTIVE

INSTRUCTIONS

Calcium

Model WAK-Ca

PC color comparison Method


Main reagent: Phthalein Complexone

Range: 0 - 50 mg Ca/L (ppm)

0 - >125 mg CaCO₃ /L (ppm)

Harmful

How to use




(1) Remove the line to clear the aperture from the top of the tube.



(2) Press the sides of the tube to expel approximately half of the volume. Maintain pressed.



(3) Immerse the tube in the sample. Release the sides to fill the tube up to the half. Shake the tube lightly a few times.



(4) After 2 minutes, put the tube on the color chart as shown and compare with the standard colors.

How to read the test

After the reaction time, compare the color of the tube with the standard colors. The nearest color indicates the measured value of the sample. A color between two standard colors indicates a value between the two standard values. Calcium hardness value is printed below the standard colors.

Care in handling of PACKTEST before and after use

Keep PACKTEST out of the reach of children.

Keep PACKTEST in a cool, dry and dark place.

PACKTEST should be thrown with burnable garbage. Conform to the legislation of waste management.

Use a package as soon as possible after opening.

The PACKTEST tube must not be opened before and after use.

Reagent is deteriorated by a temperature higher than 50°C

First Aid Measures

Contains Sodium Tetraborate. It is a harmful.

Eye contact → Immediately rinse eyes with water for at least 15 minutes. Consult a physician.

Skin contact → Immediately flush skin with water.

Ingestion → Immediately rinse mouth. Consult a physician.

In case of doubt, consult a physician.

**KYORITSU CHEMICAL-CHECK Lab., Corp.**37-11, DEN-ENCHOFU 5 CHOME, OHTA-KU, TOKYO 145-0071 JAPAN
FAX: 81-3-3721-0666 <http://kyoritsu-lab.co.jp>

PACKTEST Calcium

Cautions

1. The Calcium PACKTEST can measure only calcium ion Ca^{2+} concentration. To measure total calcium, sample must be treated to dissolve precipitated forms before measurement.
2. The calcium hardness read lower values of standard colors, or the calcium value times 2.5.
3. The normal pH range is 6 -11. If necessary, adjust the pH with diluted sulfuric acid or sodium hydroxide solution.
4. Ensure that PACKTEST tube is filled up to the half.
5. Partially undissolved reagent will not affect the measurement.
6. Keep sample temperature in the range 15°C - 40°C . Lower temperature necessitates longer reaction time.
7. Read the test under a daylight type lamp.
8. Put the line back into the aperture after using to prevent reagent spilt.

Interferences

Standard colors were determined from standard solutions. However, coexisting substances will cause inaccurate results. The list below reports substances concentrations under which ones interferences are insignificant:

$\leq 1000 \text{ mg/L}$: Cl^- , K^+ , Na^+ , NO_2^- , NO_3^- , SO_4^{2-} , Surfactant

$\leq 500 \text{ mg/L}$: PO_4^{3-}

$\leq 200 \text{ mg/L}$: NH_4^+

$\leq 20 \text{ mg/L}$: Zn^{2+}

$\leq 10 \text{ mg/L}$: Cu^{2+} , Mn^{2+}

$\leq 8 \text{ mg/L}$: Mg^{2+} ,

$\leq 1 \text{ mg/L}$: Residual chlorine

$\leq 0.5 \text{ mg/L}$: Fe^{2+} , Fe^{3+}